

**JORDAN** | RECLINERS



**krug**

HEALTHCARE SOLUTIONS



## JORDAN SLEEP RECLINER







**JORDAN 3 POSITION RECLINER**



## JORDAN BARIATRIC RECLINER







**JORDAN INFUSION RECLINER**









**JORDAN** MOTORIZED RECLINER

## OPTIONAL FEATURES



### PUSH BAR

This is a standard feature on the Jordan recliner, ergonomically designed to assist caregiving staff in the movement and operation of the recliner, and reduce back strain on the caregiver.



### CUSHIONED ARMREST

The Infusion arm was designed to maximize patient comfort during infusion treatments. The extra padded and wider arms provide comfort and support over long periods of time. Arms can be ordered in Krug standard healthcare rated gray vinyl or in coordinating patterns.



### FOLEY HOOK

This recessed, retractable hook system provides a mounting system for catheter bags. It is spring-loaded to automatically retract when not in use, to eliminate any hazard.



### IV POLE

The Jordan IV pole is height adjustable, and can be added in the field. It can be located in any of four locations (left or right, front or back) - if ordered with the recliner, the location must be specified at time of order.



### CENTRAL CASTER LOCKING

Allows for single action locking of all 4 casters at one time. Casters are locked both directionally and rotationally, keeping the recliner in place when locked. Color coded external foot pedal - Red for lock and Green for unlocked.



### COMBINATION UPHOLSTERY

Combination upholstery is an option on all Jordan Sleep Recliners.



### SWING AWAY TABLET

The tablet has the ability to slide up and fold over and rotates 180°. The break away design allows for the tablet to be pushed away from the patient quickly in case of an emergency.



### FLIP UP TABLET

The tablet is located opposite side of the ottoman lever. It is available in Beech and Maple veneer, as well as plastic laminate with a PVC edge. It folds flat (protruding 3" from the side of the recliner when flat).



### THERMAL COMFORT

Providing warmth to patients under going chemotherapy and other treatments can be vitally important. Thermal Comfort is operated with a hand-held controller located on the inside arm of the recliner. It provides two heat zones: one in the seat and one in the back. Back and seat zones can be controlled independently, at a high or low heat setting. Thermal Comfort features a safety time-out that automatically shuts off heat after 30 minutes of use.



### POWER CONTROL

Optional motorized recliner is available on Jordan Recliner and Jordan Sleep Recliner. The standard recliner features a single function control which operates the back and ottoman simultaneously. Sleep Recliners feature independent back and ottoman control.



### FOOTREST LEVER

The Footrest Lever is standard on Jordan Sleep Recliners. It extends and closes the built-in foot rest.



### PATIENT TRANSFER ARM

The Patient Transfer Arm is designed to fold down where the arm is level with the chair seat. This design allows the caregiver to easily transport patients from the chair to bed. The arm is folded down by accessing the adjustment lever at the side of the chair. Please note, the Patient Transfer Arm cannot be located on the same side as IV Poles, and tablet option or the footrest lever.



## STANDARD FEATURES



### CLEANING & INFECTION CONTROL

Specifically designed to ease the cleaning and disinfection of the recliner and the floor underneath. Upholstery covers can be made removable, and individual components are automatically removable, to enable thorough cleaning and disinfection.



### CASTER

The standard caster is 4" in diameter, soft wheel, full swivel, with a lock that in a single action locks both the swivel and rotation of the caster. The Central Locking caster has the same features, with a hooded design.



### ENHANCED CUSHIONING

The Jordan Recliners have been upgraded with a softer seat, a more pronounced headrest, and contoured back shaping for improved support. These enhancements, along with the addition of dual-density foam, make these recliners the leader in patient comfort.



### REMOVABLE COVERS AND REPLACEMENT COMPONENTS

Removable seat and back covers can be quickly installed on Jordan products that were originally ordered with the Removable Cover Option. Removable seat and back covers can also be installed on Jordan products that were not originally ordered with the Removable cover option, with some simple staple removal on the original seats and backs.



### HIGH DURABILITY MECHANISM

Jordan Recliner mechanisms are designed for high-traffic, intensive use in care-giving environments. They are built with high-strength steel, and include Delrin™ wear pads between moving components for added smoothness and longevity.



### LEG



All products are air emissions certified and are available as FSC® certified

# krug

HEALTHCARE SOLUTIONS