Jordan Compact Recliner



krug healthcare §

Jordan Compact Recliner

The Jordan Compact Recliner is lighter-scaled, lower weight and dimensionally smaller than most sleep recliners, but maintains comfortable proportions for extended sitting, resting and treatment. It is a product that has some unique solutions to very real and practical needs and the perfect solution for patient rooms and treatment areas where available space and freedom of movement are important.





Fully Upholstered Arm

Jordan Compact Recliner Arm Styles



Fully Upholstered Arm with Urethane Arm Cap



Fully Upholstered Arm with Solid Surface Arm Cap



Jordan Compact Recliner

Enhanced Cushioning

The Jordan Compact Recliners have been designed with a softer seat, a more pronounced headrest, and contoured back shaping for improved support. These enhancements, along with the addition of dual-density foam, make these recliners the leader in patient comfort.



.

The Jordan Compact Recliner features a Infinitely Adjustable Mechanism:

Position A - seated with feet on the floor, footrest in down (closed) position. Seat-back angle is 102 degrees.

Position B - footrest is fully extended. The seatback angle remains at 102 degrees.

Position C - is the range of infinite adjustability from the seated position of 102 degrees back angle to the fully flat position at a maximum of 175 degrees. The Infinitely Adjustable Back angle can be locked in any position with the back adjustment lever.

Jordan Compact Recliner Infinitely Adjustable



an

HIGH DURABILITY MECHANISM

Compact Recliner mechanisms are This mechanism provides the additional The Jordan Compact Recliner The footrest lever extends and closes designed for high-traffic, intensive use in Trendelenberg position: full recline with features a Black back adjustment lever the built-in footrest. The footrest lever is care-giving environments. They are built footrest fully extended, the back is now which creates an infinite number of finished in powder-coated Silver to with high-strength steel, and include 65 degrees from vertical - resulting in back adjustments within the chair. compliment the casters. Delrin™ wear pads between moving the feet being above the head. components for added smoothness and longevity.



TRENDELENBERG MECHANISM



BACK ADJUSTMENT LEVER

Conveniently located on the care-giver side of the arm.



FOOTREST LEVER



FOOT TRAY

when standing on the tray.



PATIENT TRANSFER ARM

down by accessing the adjustment cycles on a full charge. lever at the side of the chair.



POWER CONTROL

This option provides a concealed, pull-out The Patient Transfer Arm is designed Optional motorized recline features a Providing warmth to patients undergoing tray to assist patients with ingress and to fold down where the arm is level with single function control which operates chemotherapy and other treatments egress. It features a non-slip tread and the chair seat. This design allows the back and ottoman control independently. can be vitally important. Thermal Comfort dedicated casters that prevent tipping caregiver to easily transport patients. Power also can also be specified with a features a recessed controller which from the chair to bed. The arm is folded lithium battery pack able to run 200 provides two heat zones: one in the



THERMAL COMFORT OPTION

seat and one in the back, including a safety time-out feature after 30 minutes.



CASTERS

soft wheel, full swivel, featuring a lock locking of all four casters at one time. assist care-giving staff in the movement that in a single action locks both the Casters are locked both directionally and operation of the recliner, and portion of the upholstered ottoman, swivel and rotation of the caster.



CENTRAL CASTER LOCKING

The standard caster is 4" in diameter, This option allows for single action Pushbar is ergonomically designed to and rotationally, keeping the recliner in reduce back strain on the care-giver. place when locked.



PUSH BAR



FOOTREST PROTECTIVE COVER

This optional clear protective cover is sewn into (and over top of) the lower providing a layer of protection. It features an anti-fog additive to maintain clarity, and is resistant to cleaning solutions.



SWING AWAY TABLET

quick access to the patient.



FLIP UP TABLET

from the user in an emergency to gain Maple veneer, as well as Laminate with a the recliner. Polymer edge.



IV POLE & HOLDER



FOLEY HOOK

The Swing Away Tablet is designed to The Flip up tablet folds flat (protruding 3" The Jordan IV Pole is height adjustable This recessed, retractable hook system be easily operated by the caregiver or from the side of the recliner when flat) (in a range of 27"), and can be added in provides a mounting system for catheter user. The Break Away design allows for and it can be located on either side of the field. It can be located on either the bags. It is spring-loaded to automatically the tablet to be easily pushed away the recliner. It is available in Beech and left or right side (facing) at the back of retract when not in use, to eliminate any hazard.

STATEMENT OF LINE



JOR8-CSRUPHS1







JOR8CSRUURS1





JOR8-PCSRUURS1





JOR8-CSRTUURS1



JOR8-PCSRTUPHS1



JOR8-PCSRTUURS1

URETHANE COLOR OPTIONS







GROUP GRADE A - COLOR

Taupe

SOLID SURFACE

GROUP GRADE A - WHITE

Cameo White





Bisque

Dove





Glacier White

GROUP GRADE B









Deep Night Sky



Jordan Compact Recliner Environmental

The use of recycled and recyclable materials, combined with environmentally-friendly production processes and packaging, have made Faeron a levelTM-certified product, and a contributor to LEEDTM-driven specifications. Faeron is certified by the GreenguardTM standard for indoor air emissions, and it's unique capability for extending product life makes it a responsible environmental choice - and one that will provide exceptional return on investment over time.







All products are air emissions certified and are available as FSC® certified

