## Jordan Compact Recliner



# **krug** HealthCare **\***

## Jordan Compact Recliner

The Jordan Compact Recliner is a lightly-scaled version of our award-winning Jordan Sleep Recliner. Dimensionally smaller than most sleep recliners, the Jordan Compact Recliner maintains comfortable proportions for extended sitting, resting and sleeping. It is an outstanding solution for areas where available space and freedom of movement are important.





Fully Upholstered Arm



Fully Upholstered Arm with Urethane Arm Cap

Jordan Compact Recliner Arm Styles



Fully Upholstered Arm with Solid Surface Arm Cap





## Jordan Compact Recliner Enhanced Cushioning

The Jordan Compact Recliner is designed with a contoured back, pronounced integrated headrest and spring seat construction supplemented with dual density foam for superior support. These enhancements make this recliner a leader in patient comfort.



## Jordan Compact Recliner Versatile Functionality



Upright Seated

Extended Footrest

## Infinitely Adjustable



## HIGH DURABILITY MECHANISM

Delrin<sup>™</sup> wear pads between moving the feet being above the head. components for added smoothness and longevity.



### TRENDELENBURG MECHANISM

Compact Recliner mechanisms are This mechanism provides the additional The Jordan Compact Recliner The footrest lever extends and closes designed for high-traffic, intensive use in Trendelenburg position: full recline with features a Black back adjustment lever the built-in footrest. The footrest lever is care-giving environments. They are built footrest fully extended, the back is now which creates an infinite number of finished in powder-coated Silver to with high-strength steel, and include 65 degrees from vertical - resulting in back adjustments within the chair. compliment the casters.



### BACK ADJUSTMENT LEVER

Conveniently located on the care-giver side of the arm.



### FOOTREST LEVER



### FOOT TRAY

This option provides a concealed, pull-out The Patient Transfer Arm is designed Optional motorized recline features a Providing warmth to patients undergoing tray to assist patients with ingress and to fold down where the arm is level with single function control which operates chemotherapy and other treatments egress. It features a non-slip tread and the chair seat. This design allows the back and ottoman control independently. can be vitally important. Thermal Comfort dedicated casters that prevent tipping caregiver to easily transport patients. Power also can also be specified with a features a recessed controller which when standing on the tray. from the chair to bed. The arm is folded lithium battery pack able to run 200 provides two heat zones: one in the down by accessing the adjustment cycles on a full charge. lever at the side of the chair.



### CASTERS

soft wheel, full swivel, featuring a lock locking of all four casters at one time. assist care-giving staff in the movement that in a single action locks both the Casters are locked both directionally and operation of the recliner, and portion of the upholstered ottoman, swivel and rotation of the caster.



### CENTRAL CASTER LOCKING

The standard caster is 4" in diameter, This option allows for single action Push bar is ergonomically designed to and rotationally, keeping the recliner in reduce back strain on the care-giver. place when locked.



### PUSH BAR



### FOOTREST PROTECTIVE COVER This optional clear protective cover is

sewn into (and over top of) the lower providing a layer of protection. It features an anti-fog additive to maintain clarity, and is resistant to cleaning solutions.



# SWING AWAY TABLET quick access to the patient.



### PATIENT TRANSFER ARM



### POWER CONTROL



### THERMAL COMFORT OPTION

seat and one in the back, including a safety time-out feature after 30 minutes.



### FLIP UP TABLET

The Swing Away Tablet is designed to The Flip up tablet folds flat (protruding 3" The Jordan IV Pole is height adjustable This recessed, retractable hook system from the user in an emergency to gain Maple veneer, as well as Laminate with a the recliner. Polymer edge.



### IV POLE & HOLDER



### FOLEY HOOK

be easily operated by the caregiver or from the side of the recliner when flat) (in a range of 27"), and can be added in provides a mounting system for catheter user. The Break Away design allows for and it can be located on either side of the field. It can be located on either the bags. It is spring-loaded to automatically the tablet to be easily pushed away the recliner. It is available in Beech and left or right side (facing) at the back of retract when not in use, to eliminate any hazard.

### 10

### STATEMENT OF LINE



JOR8-CSRUPHS1



JOR8CSRUURS1



JOR8-PCSRUPHS1












JOR8-PCSRUURS1

JOR8-CSRTUURS1 JOR8-CSRTUPHS1

JOR8-PCSRTUPHS1

### URETHANE COLOR OPTIONS



GROUP GRADE A - COLOR

### SOLID SURFACE

GROUP GRADE A - WHITE



GROUP GRADE B







GROUP GRADE C

Bisque







Antarctica

Savannah

Dove

Deep Night Sky





## Jordan Compact Recliner Environmental

The use of recycled and recyclable materials, combined with environmentally-friendly production processes and packaging, have made Jordan a level<sup>™</sup>-certified product, and a contributor to LEED<sup>™</sup>-driven specifications. Jordan is certified by the Greenguard<sup>™</sup> standard for indoor air emissions, and it's unique capability for extending product life makes it a responsible environmental choice - and one that will provide exceptional return on investment over time.



All products are air emissions certified and are available as FSC® certified



solutions@krug.ca | www.krug.ca | 1.888.578.KRUG